

Pharmacy Refrigerator

HYC-610



LED Temperature Display

Display precision 0.1°C.



Multiple Alarm Functions

Alarms including high and low temperature, remote alarm, sensor error, power failure, low battery, door ajar.



Adjustable Shelves

The shelves are adjustable to any height to fulfill different user requirements.



Large Storage Space



Lockable Design



Self Closing Door with 90° Stay Open Feature

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International



Product Advantages



Energy Efficient Hydrocarbon (HC) Refrigeration System

Imported efficient inverter compressor with optimized cooling system, more energy efficient.



Greener by Design

Hydrocarbon compressor and refrigerant, the refrigerant does not contain any fluorine or chlorine, greener and environmentally friendly.



Energy conserving, maintains temperature uniformity to protect valuable samples.



Superior cabinet temperature uniformity, automatic defrosting



Specifications

Model	HYC-610				
Technical Data	Cabinet Type		Upright, Single Door		
	Climate Class		N		
	Cooling Type		Forced Air Cooling		
	Defrost Mode		Auto		
	Refrigerant		HC		
	Sound level (dB (A))		43	45	
Performance	Temperature Range (°C)		2~8		
Control	Controller		Microprocessor		
	Display		LED		
Electrical Data	Power Supply (V/Hz)		220~240/50/60		115/60
	Electrical Current (A)		3	6.2	
Dimensions	Capacity (L/Cu.Ft)		610/21.6		
	Net/Gross Weight (approx)	kg	204/227		
		lbs	449.7/500.4		
	Interior Dimensions (W*D*H)	mm	680*640*1400		
		in	26.8*25.2*55.1		
	Exterior Dimensions (W*D*H)	mm	780*840*1945		
		in	30.7*33.1*76.6		
	Packing Dimensions(W*D*H)	mm	885*965*2110		
		in	34.8*38.0*83.1		
Container load (20'/40'/40'H)		12/26/26			
Alarms	High/Low Temperature		Y		
	Remote Alarm		Y		
	Power Failure		Y		
	Sensor Error		Y		
	Low Battery		Y		
	Door Ajar		Y		
Accessories	Caster		Y		
	Foot		Y		
	Porthole		Y		
	Shelves/Drawers		6/-		
	USB Interface		Optional		
	Printer		/		
	Temperature Recorder		Optional		
Others	Certification		CE	UL	

*Haier Biomedical reserves the right to change products and specifications without prior notice.